


Tru-Vu[®]


Mini Hysteroscope



5729B - 2.9 mm x 30 Degree
Tru-Vu[®] Rod Lens Hysteroscope
Steam Autoclavable
Color Code: **Red**

2600B - 4.3 mm Inner Sheath and Bridge,
5 FR Channel for Operating Instruments
Continuous Flow

2600S - 5.0 mm Outer Sheath
with Continuous Flow Valve
Tear Drop Shape for
Easy Passage



FEATURES:

- Small Diameter
- Continuous Flow System
- 5 FR Instrument Channel
- Tru-Vu[®] Precision Rod Lens Endoscope
- High Aperture Fiberoptics Light Transmission
- Compatible with *Essure[™] Procedure
- Tear Drop Design for Easy Passage
- 100 % Stainless Steel Construction
- Ease of Reprocessing - Cold Soak - ETO Gas - Steam Autoclavable



2600C - 5 FR Blunt Scissor
Semi Rigid - Single Action



2600P - 5 FR Punch Forcep
Thru Cutting - Semi Rigid



2600K - 5 FR Biopsy Spoon Forcep
Semi Rigid

*Essure[™] is a Conceptus (CPTS) Procedure



**Advanced
Endoscopy
Devices**

Endoscopy Specialists since 1985

22134 Sherman Way
Canoga Park, CA 91303 - U.S.A.
phone: 818-227-2720 - fax: 818-227-2724
www.aed.md
E-mail: advendodev@earthlink.net