



Advanced
Endoscopy
Devices

Precise Lap Bulldog Clamp Forcep



Model # PB500

- Lap Forcep for Applying Clamp
- 12.5mm x 38cm
- Rotatable Shaft
- Cleaning Flush Port on Rotator
- Spring Action Handle
- Included with a 45mm and 60mm Clamp



PB545

45mm Clamp



PB560

60mm Clamp

1-818-227-2720

Advanced Endoscopy Devices, Inc.
22134 Sherman Way, Canoga Park, CA 91303 USA
www.aed.md sales@aed.md